

Remote Monitoring for Medtronic Heart Device Patients

Sending device information
to your clinic without
an analog telephone landline
just got simpler.

If you want to use the CareLink Patient Monitor to send information stored in your heart device to your doctor, but do not have access to an analog telephone landline, we may have a solution for you.

You can now send your information using the WireX™ cellular accessory.

The accessory connects to your CareLink Patient Monitor and uses cellular signals to send your device information to your clinic. It's a simple way to stay connected to your doctor via the CareLink Network.

NOTE: Lack of cellular data coverage is one of the major factors that can contribute to unsuccessful use of the WireX™ cellular accessory. Environmental factors such as cell tower saturation, landscape topography, home construction materials, and in-home electrical interference may result in connectivity issues that prevent the use of cellular signals to transmit data to the CareLink Network. Based on the nature of cellular data communication, the WireX accessory may not be a solution for some patients.

**For more information on how to order
and use the WireX cellular accessory,
call 1 (877) 609-6698.**

Additional Device Information

The Medtronic CareLink® Monitor is a prescription device indicated for use in the transfer of patient data from some Medtronic implantable cardiac devices based on physician instructions and as described in the product manual. This product is not a substitute for appropriate medical attention in the event of an emergency and should only be used as directed by a physician.

The CareLink® Service is prescribed by your physician. This service is not for everyone. Please talk to your doctor to see if it is right for you. Your physician should discuss all potential benefits and risks with you. Although many patients benefit from the use of this service, results may vary. For further information, please call the Medtronic toll-free number at 1 (800) 929-4043 (Monday-Friday, 8:00 a.m. to 5:00 p.m., Central time) or see the Medtronic website at www.medtronic.com.

www.medtronic.com

World Headquarters
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel: (763) 514-4000
Fax: (763) 514-4879

Medtronic USA, Inc.
Toll-free: 1 (800) 328-2518
(24-hour technical support
for physicians and medical
professionals)



Connections. Simplified.

CELLULAR ACCESS TO THE
MEDTRONIC CARELINK® NETWORK



UC20090505c4 EN © Medtronic, Inc. 2014. Minneapolis, MN. All Rights Reserved. Printed in USA. 01/2014

Benefits of the WireX™ cellular accessory:

Convenience

The WireX cellular accessory allows you to send your device information to your clinic, as instructed by your doctor, without an analog telephone landline or cell phone.

Affordability

Compared to costs associated with in-person doctor appointments, digital and analog telephone landline charges, this is an affordable solution for remote monitoring of your heart device.*

Mobility

The cellular accessory offers you freedom to travel globally and stay connected to your doctor via the CareLink® Network.**

The CareLink Network: Now offering another way for you to stay connected to your clinic – anytime and anywhere – with the WireX cellular accessory.

* Medtronic data on file.

** Please call 1 (877) 609-6698 to confirm cellular signal coverage for a specific location.

Common Questions and Answers

What is the WireX cellular accessory?

The WireX cellular accessory connects to the CareLink Patient Monitor and uses cellular signals to send information stored in your heart device to your clinic via the CareLink Network. **It will only work with the CareLink Patient Monitor.**

How do I know if I need the cellular accessory?

If you have a phone service that is not compatible with the CareLink Patient Monitor, you do not have an analog telephone, or you frequently travel globally, you may benefit from the cellular accessory.

Can I use the cellular accessory when traveling?

Yes, you can use the accessory with your CareLink Patient Monitor to send device information to your clinic when traveling within or outside of the United States. However, cellular signal coverage may not be available in all locations. You can call Medtronic at 1 (877) 609-6698 to find out if coverage is available in a certain location.

Do I need a cellular phone to use this accessory?

No, you do not need a cellular phone to use the accessory.

How much does the cellular accessory cost?

You will be charged a monthly service fee of \$10.99 per month, plus any applicable taxes, regardless of how often you use it. Medtronic pays for the cellular accessory hardware.

How do I order the cellular accessory?

After your clinic has ordered a CareLink Patient Monitor for you, you will need to call Medtronic at 1 (877) 609-6698 to order the cellular accessory. Further information on payment and cellular signal coverage will be discussed with you when you call to place your order.


Can I use my existing cellular phone plan to send device information to my clinic?

No, Medtronic only offers the WireX cellular accessory as a reliable and secure way for you to send device information to your clinic using cellular signals.

How does the WireX cellular accessory connect to the CareLink Patient Monitor?

The accessory comes with a pre-connected phone cord that plugs into your CareLink Patient Monitor. The accessory plugs into a wall outlet for power.

Battery-Operated CareLink Patient Monitor



Electricity-Powered CareLink Patient Monitor

